

# molecules to man imaging

## bioVet

Our market leading accessory portfolio of high quality exceptionally performing products will give you the workflow solutions you need to drive your research to higher levels. We're redefining the gold standard in utilization of accessories to help you produce better results, more reliably, with higher throughput: faster.

What is your requirement? Our accessories are designed with researchers to optimize workflow and dramatically reduce the time you spend in pre and post experiment preparation. Talk to us to learn why our preclinical multi modal accessories are sold by the world's largest brands in the industry: Siemens Molecular Imaging, GE Healthcare, Philips Medical Systems, GammaMedica Ideas, Molecular Imaging Products Company and Scanco to name a few.



**m2m's BioVet** comprehensive and fully featured system includes everything you need to get started, and then develop and extend your workflow and scanning protocols.

### About Us:

**Join the trend of researchers choosing us**

*YOU are our luminary customer base!*

Achieve otherwise unobtainable image quality and drive your research to new levels of discovery with better products, workflow, and customer support. Our cutting edge products are designed, developed and manufactured to strict ISO quality standards.

Learn more about why more and more researchers are turning to us to improve their results and workflow: Talk to us today about how we can help you.

### BioVet Physiological Monitoring and Gating:

**Gold standard in multi modal gating**

- Our preclinical monitoring devices are the market leading solution in PET, CT and SPECT preclinical imaging. BioVet's monitoring capabilities and automation features allow advanced triggering schemes for use with market leading instruments.
- Monitors: ECG, Respiration and Temperature
- Optional subject heating channel
- USB powered with dual and single trigger outputs
- Easy to use: software designed to make the setup process fast and accurate
- Online support and product updates (firmware and software)



## bioVet accessories

### **Multi Modal Animal Beds: Standardized platforms for preclinical imaging**

Developed in conjunction with UCLA's Crump Molecular Imaging group, our multi modal animal beds are optimized for co-registration between imaging modalities and high throughput. We have solutions available for many platforms including Siemens Inveon, Concorde and Imtek systems (now all Siemens products).

- Integrated subject heating
- Standard platform between systems
- Anesthesia delivery and scavenge
- Mounting plates to allow translation between systems without removing the subject
- Sterilization compatible

### **Integrated and Stand Alone Heating Solutions:**

Heating is well known to be an important to keep the subject comfortable — and to ensure your consistent metabolic results. Our solutions utilize an imaging compatible resistive heating mat which is available in various sizes.

- Heating solutions to work with the BioVet system
- Stand alone heating solutions utilizing a small controller with LCD screen to control temperature
- Heater mats: different sizes for different applications
- No hot air and no warm water — No leakage and no dehydration

### **Global coverage and service: We're bringing clinical levels of service to preclinical markets**

We offer full line service and extended warranty on all our products. Speak to us today to find out how to minimize downtime through service agreements with us. With offices in North America and Asia-Pacific, you can be confident that someone is always available to support you.



190 Alpha Park  
Cleveland, OH 44143 USA  
+1 877 887 8989 (toll free)  
+1 440 684 9695 (fax)  
sales@m2mimaging.com  
www.m2mimaging.com