

molecules to man imaging

BioVet Multimodal Chambers

Our market leading accessory portfolio is setting the standard for high performance, cost effective, and workflow enhancing solutions for academic, pharma and drug discovery, and government labs globally. Our products are developed in close collaboration with researchers, technicians, and equipment manufacturers from around the world.

What is your requirement? Our accessories are designed with your labs productivity and efficiency in mind. Careful attention is paid to workflow, pre and post experiment setup and close down, and longitudinal reproducibility of experiments over the long term.



m2m's BioVet Multimodal Chambers are a comprehensive and fully featured system which includes everything you need to get started, and then grow with you as you develop applications and workflow.

About Us:

Join an increasing number of researchers and manufacturers who choose our products

You are our luminary customer base!
Achieve otherwise unobtainable image quality and drive your research to new levels of discovery with better products, improved workflow, and world class support. Our products are designed and manufactured to ISO quality standards.

Mouse Imaging System includes:

- The Chamber, with fiducial pockets
- Head First Neuro Sled with sliding nosecone and bitebar
- Tail First Sled & Syringe Holder
- Mounting adapter
- Heater pad
- Manual

Multi Modal Imaging Chambers: Standardized platforms for preclinical imaging

Our imaging systems are designed as detachable chambers with application specific insert sleds. This allows for efficient setup and prep away from the imager. Our nosecones are specifically designed for optimum performance, also provide exhaust management, and positive closure.

Materials and shapes are carefully chosen for minimal attenuation, no artifacts, and ease of cleaning.

All our chambers are designed for use with our BioVet physiological monitoring and warming products.

Developed in conjunction with UCLA's Crump Molecular Imaging Group, they are optimized for positioning control and reproducibility, co-registraion between modalities, and provide application specific positioning for optimum image quality.

Global coverage and service

We're bringing clinical levels of service to preclinical markets

We offer full line service and warranty support on all our products that you can rely on. Just ask our customers!



BioVet

Multimodal Chambers: Mouse Imaging System



Tail First Mouse Sled



Head First / Neuro Mouse Sled

