

molecules to man imaging

microimaging

Our market-leading Microimaging product portfolio of high value products will give you more tools per budget dollar than any other supplier. Integrated systems, probe bodies and detachable probes, in vivo imaging platforms, our highly uniform B1 field SAW platform, and single, double, triple tuned, and broadband designs, we can create a solution for your specific research needs on any OEM platform.



What is *your* requirement?

Our product portfolio consists of thousands of designs ranging from 2mm to 30cm, from 100 to 900 MHz, every nmr visible nuclei, and for Bruker and Varian platforms. Chances are we've built it before and if not, we'll build it for you...

Our solutions for any field strength include:

- Probes and probe bases
- volume coils
- surface coils
- ex vivo sample probes
- integrated in vivo platforms
- Micro 2.5 imaging probes, MicroMouse 2.5, Micro 2.5 AHS/RF Probes
- Mini 0.36 systems, Mini 0.36 AHS/RF
- Mini 0.5 systems, Mini 0.5 AHS/RF
- Widebore and Super Widebore probes

▶ **Our probe bases have clutched cap drives. Never break a variable cap again.**

About Us

All our customers are luminary users to us!

Achieve otherwise unobtainable image quality while getting the most for your research budget expense. Use our high quality reliable products as the daily tools of your research. Our workflow enhancing platform designs let you get the most scanning done in a day of any products available. Our products are developed and manufactured to global ISO 9001:2008 standards.



microimaging

Product Portfolio

Choose from any of our over 5000 designs of microimaging probes and accessories, or set a specific requirement and we will develop a unique solution just for your project. We have been providing serial products and platform systems, end user custom products, and private label solutions for OEM's for decades.

Accessories and Interfacing

We make a wide variety of active and passive combiners, and universal platform agnostic coil control units, "CCU's" that function with all OEM platforms. We can provide products with the new 'monoplug' clinical style connectors as well.

Repairs, Service, and Upgrades or Modifications

We service what we build and sell, (and we can also service what others do...). In 90% of cases, our repair turnaround times are 48 working hours or less. Don't be without your valuable research tools for weeks or months at a time while they disappear without notice into the repair cycles of the system manufacturers. We can also frequently take existing probes and provide upgrades or modifications if your research needs have changed. Contact us for estimates.



+1 440-684-9690 (office)
+1 440-684-9695 (fax)
+1 877-887-8989 (toll free)

sales@m2mimaging.com
info@m2mimaging.com
support@m2mimaging.com
www.m2mimaging.com